





simpas

PTCA BALLOON DILATATION CATHETER


P L U



No Lesion Can Stop You
Unique Lesion Crossability
with

■ Excellent Hydrophilic coating

■ Low Balloon Crossing Profile



- RX PTFE Coated Hypotube Design
- Excellent Hydrophilic coating
- 0,016" Tip Profile
- Enhanced Pushability
- Enhanced Average
- Atraumatic Tapered

Balloon Diameter (mm)	Balloon Length (mm)					
	10	12	15	20	25	30
1.25*	PS12510-T	---	PS12515-T	PS12520-T	–	–
1.50*	PS15010 PS15010-T	PS15012 PS15012-T	PS15015 PS15015-T	PS15020 PS15020-T	–	–
2.00*	PS20010 PS20010-T	PS20012 PS20012-T	PS20015 PS20015-T	PS20020 PS20020-T	PS20025	–
2.25	PS22510	PS22512	PS22515	PS22520	PS22525	PS22530
2.50	PS25010	PS25012	PS25015	PS25020	PS25025	PS25030
2.75	PS27510	PS27512	PS27515	PS27520	PS27525	PS27530
3.00	PS30010	PS30012	PS30015	PS30020	PS30025	PS30030
3.50	PS35010	PS35012	PS35015	PS35020	PS35025	PS35030
4.00	PS40010	PS40012	PS40015	PS40020	PS40025	PS40030
4.50	PS45010	–	PS45015	PS45020	PS45025	PS45030
5.00	PS50010	–	PS50015	PS50020	–	–
6.00	PS60010	–	PS60015	PS60020	–	–

* Single Marker Band is available for 1.5 mm & 2.0 mm diameters.

